Quick protein binding analysis by label-free thermal shift analysis on the Tycho NT.6

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
Figure 1: Monitoring p38 kinase binding activity using a label-free thermal shift analysis. A) Interaction of 4 μ M p38 α kinase with 5 mM Mg²⁺ ions. B) Interaction of 4 μ M p38 α kinase with 500 nM SB203580. Vertical lines indicate inflection temperatures (T_i), colored circles on the y-axis indicate initial ratio values.


Figure 2: Label-free thermal shift analysis of 10 μM maltose-binding protein (MBP) interacting with 100 μM maltose using Tycho NT.6. Dramatic differences in the initial ratio values are likely due to the conformational changes of MBP upon binding of maltose. Vertical lines indicate inflection temperatures (T_i), colored circles at the y-axis indicate initial ratio values.


Figure 3: Testing the interaction of bacterial cell wall complex (NAG) with lysozyme using the Tycho NT.6. 5 μM lysozyme was incubated with 50 μM 3'NAG and analyzed. Vertical lines indicate inflection temperatures (T_i), colored circles at the y-axis indicate initial ratio values.

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Reference

1. Telmer, P.G. and B.H. Shilton, Insights into the conformational equilibria of maltose-binding protein by analysis of high affinity mutants. *J Biol Chem*, 278(36): p. 34555-67 (2003).

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NanoTemper Technologies GmbH, Floessergasse 4, 81369 Munich