



Proteus Ni Advance Spin Column Kit User Guide

Protein Ark's HiFLiQ Ni Advance spin column kit, is designed for rapid one-step batch binding purification, and ideal for preparative purification and contaminant removal. The Ni Advance FPLC spin column kit is supplied as loose resin with mini or midi spin columns, designed for affinity purification of poly-histidine tagged recombinant proteins by immobilized metal ion affinity chromatography (IMAC).



The Ni Advance Agarose resin provides high binding capacity with minimal Ni²⁺ ion leakage for high stability, chemical compatibility and reuse. Compatible with all common bench top and micro centrifuges which accept both 50ml and 2ml tubes.

Features of the Ni Advance spin column kit:

- Fast and reliable affinity purification.
 - Nickel ions are carefully loaded onto a 6% highly cross-linked agarose matrix via chelating coupled ligand to obtain a stable affinity matrix with a high binding capacity for histidine residues.
 - High protein binding capacity combined with minimal Ni²⁺ ion leakage.
 - Keep your buffer in its preferred conditions as Fastback Ni Advance resin is resistant to EDTA (up to 20mM) and DTT (up to 20mM) for IMAC chromatography – Ni ions remain bound to ligand
 - Simple recharging procedure with Ni²⁺SO₄.
 - Rechargeable with other compatible metal ions (Co²⁺, Cu²⁺, Zn²⁺, Fe³⁺, and Al³⁺) for different binding affinities.
-

Resin Specification:

Resin:	Fastback Ni Advance
Base Matrix:	6% cross-linked Agarose
Coupled ligand:	Chelating ligand
Metal ion capacity:	>75 $\mu\text{mol/ml}$
Typical Binding Capacity:	80 mg/ml
Mean Bead Size:	90 μm
Storage (2-8°C):	20% Ethanol

Contents of Ni Advance Mini Kits:

Proteus Kit code:	NiAdv2	NiAdv 12	NiAdv 48
• Vials containing 0.25 ml Fastback Ni-Advance resin	2	12	48
• Proteus spin column barrels with clear spin push cap (20ml capacity in a swing bucket rotor)	2	12	48
• Yellow batch incubation screw caps	2	12	48
• 50 ml centrifuge tubes.	2	12	

Contents of Ni Advance Micro Kits:

Proteus Kit code:	NiAdv2Micro	NiAdv24Micro	NiAdv100Micro
• Fastback Ni-Advance resin	1 x 0.2 ml	1 x 2.4 ml	1 x 10 ml
• Proteus spin column (600 μl capacity in a fixed angle rotor)	2	24	100
• 2.2 ml centrifuge tubes	4	48	100

Chemical compatibility:

Buffer compatibility	
Phosphate buffers:	Up to 100mM phosphate buffer recommended
Tris, HEPES & MOPS:	Up to 100 mM (secondary or tertiary amines may reduce the metal ions)
pH range:	3.0-12.0
Chelating agents	
*EDTA:	Up to 20 mM
Sulfhydryl reagents	
β-mercaptoethanol:	Up to 20 mM (can cause some reduction of metal ions)
**DTT:	Up to 20 mM
TCEP:	Up to 20 mM
Denaturants	
Urea:	8 M
Guanidinium hydrochloride:	6 M
Detergents	
DM (n-Decyl-β-Dmaltopyranoside):	1.00%
DDM (n-Dodecyl-β-Dmaltoside):	1.00%
NM (n-Nonyl-β-Dglucopyranoside):	1.00%
OG (n-Octyl-β-Dglucopyranoside):	1.50%
TDM (n-Tetradecyl-β-Dmaltopyranoside):	0.005%
Triton®:	2.00%
Tween®:	2.00%
NP-40:	2.00%
Cymal 6:	1.00%
Fos-Choline 16:	0.05%
CHAPS:	up to 1%
Other additives	
Imidazole:	500 mM
NaCl:	2.0 M (recommended concentration 300 mM)
MgCl ₂ :	Up to 4 M
CaCl ₂ :	Up to 5 mM
Glycerol:	Up to 50%
Methanol:	100%
Ethanol:	100%
Acetonitrile:	30% (v/v)

*NOTE: Stability measured after 24 hours incubation with 20 mM EDTA results in no decay in binding capacity.

**NOTE: Stability measured after 24 hours incubation with 20 mM DTT results in no decay in binding capacity.

Storage Conditions:

Remove the Fastback Ni-Advance resin vials from the kit and store it at 2-8°C. There is no need to place the rest of the kit in a refrigerator or cold room. Do not freeze the resin vials or store them at room temperature. Freezing the suspension will damage the agarose beads. The resin is pre-swollen and defined. Proteus spin columns are stable for up to 2 years at 2-8 °C from the date of manufacture. The expiration date is recorded clearly on the outside of the pack. All resin is stored in 20% v/v ethanol.

Principles of metal chelate chromatography:

Immobilized metal ion affinity chromatography (IMAC) technology was introduced by Porath *et al* (1975). Certain amino acids such as histidine, tryptophan, cysteine and tyrosine can act as electron donors on the surface of the protein and bind reversibly to a transition metal ion. IMAC purification utilizes this metal ion binding by coupling a chelating group (such as Nitrilotriacetic Acid (NTA) or Iminodiacetic Acid (IDA)) to a stable chromatographic resin (such as Agarose). Transition metal ions such as Ni²⁺, Co²⁺, Cu²⁺, or Zn²⁺ (see Porath and Olin, 1983; Porath, 1988; Sulkowski, 1989) can then be loaded and immobilised onto the chelating group enabling high affinity binding, in the majority of cases this is via an 6-8x polyhistidine-tag engineered at the N- or C-terminus of a recombinant protein (K_d -10⁻¹³ at pH 8.0). Nickel and Cobalt are the most widely used metal ions as they confer the highest affinity (Ni²⁺>Co²⁺) with the highest specificity (Co²⁺>Ni²⁺) for the IMAC tags.

The simplicity of IMAC technology is extremely attractive as it lends itself to a bind-wash-elute mode of operation if the appropriate buffer formulations are selected. IMAC purification can also be achieved using samples without any prior treatment (e.g. buffer exchange steps). The use of metal chelate affinity is widespread for the selective adsorption of engineered recombinant proteins and has largely superseded non-affinity methods of chromatography for purifying recombinant proteins.

Often a researcher needs to choose the correct base matrix, select the optimal activation chemistry and couple a suitable ligand to develop an affinity column to purify their target protein. Protein Ark's HiFliQ columns remove the guesswork associated with optimizing the resin chemistry by offering high quality 1 and 5 ml FPLC columns pre-packed and pre-charged with high capacity Fastback Ni Advance resin that satisfied stringent quality control to guarantee reproducible purification performance.

Fastback Ni Advance resin:

Nickel ions are carefully loaded onto an agarose matrix via chelating coupled ligand to obtain a stable affinity matrix with a high binding capacity for histidine residues. The resin is charged with Ni²⁺ to give a marine blue appearance. Fastback Ni Advance is very stable and can resist the following conditions in most situations: Buffers at pH 4-9, 100% methanol, 100% ethanol, 8 M urea, 6 M guanidinium hydrochloride, 30% (v/v) acetonitrile, 20 mM DTT, 20 mM EDTA. Fastback Ni Advance is also stable at pH >9.0 and can be regenerated by alkaline solutions, such as sodium hydroxide.

Application drivers for Metal Chelate chromatography:

Conditioned media secreting target proteins
Screening expression clones for high levels of His-tagged proteins.

Purification of recombinant proteins for raising antibodies.

Purification of recombinant proteins for activity and/or structural studies.

General considerations for selecting optimal binding conditions for the IMAC resin:

This resin exploits the polyhistidine sequence that permits efficient purification of the expressed protein from a broad host such as: bacterial cells, yeast cells, insect cells, or mammalian cells. Lysis conditions are dependent upon the type of expression system and nature of the recombinant protein. Mammalian or Baculovirus-infected insect cells can be lysed by sonication at +4°C with either freeze/thaw cycles or the addition of up to 1 % non-ionic detergents. While, cell lysis of *E.coli* bacteria is usually achieved by sonication on ice or homogenization either with or without lysozyme treatment.

The culture pellet containing the expressed protein is harvested and resuspended in lysis buffer at a pH close to pH 7.4-8.0 using a similar concentration of buffer, imidazole, and salt to that of a pre-equilibration buffer used for the IMAC column. Typically, a protease inhibitor cocktail, such as Boehringer "Complete EDTA-free", along with other additives like DNase I (5-50 µg/ml), Benzonase® and 10 mM β-mercaptoethanol are added to the lysis buffer. The inclusion of β-mercaptoethanol to the lysis, binding, wash and elution buffers is optional and depends upon whether the His-tagged protein elutes with any additional contaminants which may be linked via disulphide bonds. It is recommended that β-mercaptoethanol be excluded during the initial purification trials.

If the binding efficiency is found to be poor and the lysis buffer differs significantly from the pre-equilibration buffer, it is recommended that the lysate be dialysed, titrated with a concentrated stock solution, or buffer exchanged using an ultrafiltration device with a more appropriate pre-equilibration buffer.

It is imperative that the lysate is completely clear prior to loading on the column as any particulate matter (e.g. cell debris) may partially foul and clog-up the resin resulting in an increased back pressure and reduced flow rates. This will significantly increase the binding, washing and elution times and effect the final purify of the eluted protein. It is recommended that the cleared lysate be filter just prior to loading even if it has been previously filtered several days before. Ideally samples should be processed rapidly and, if the need arises, the protein purified at 4°C. It is also recommended that number of freeze/thaw cycles be minimised during storage to reduce the amount of aggregation/precipitation of the proteins.


Optimal buffer conditions for binding the target molecule to a resin are critical for successful separation of the protein. If the binding conditions are not optimal with respect to pH, salt concentration, presence of metal ions or chelating agents, flow rates, residence time ...etc, purification can be adversely affected.

Choosing the correct buffer conditions for IMAC separations:

We recommend sodium phosphate buffers as buffers with secondary or tertiary amines e.g. Tris buffers can significantly reduce the Nickel and Cobalt ions. pH levels of 7-8 works well for most immobilized Ni²⁺ applications and 0.15-0.5 M NaCl is often added to the buffers to prevent non-specific ionic interactions occurring and to stabilize the proteins in solution.

When a recombinant protein is expressed at high levels in *E.coli*, the protein can often elute as insoluble aggregates called inclusion bodies. Non-ionic detergents such as 8 M urea and 6 M guanidinium hydrochloride do not interfere with metal chelate affinity separations. These denaturants completely unfold the target protein making the His-tag much more accessible for interaction with the immobilized Ni²⁺ matrix (see the section 'protein purification under denaturing conditions' for further details).

Batch Binding Protocol:



Mini Kit Protocol

Recommended Protocol:

The following spin speeds and times are appropriate for a 0.25–1ml resin bed volume. Spin times for each of the following steps may increase with larger bed volumes.

PRE-EQUILIBRATION

1. Remove the CLEAR spin push cap and pipette 0.5 ml Ni-Advance resin slurry (50% slurry: 0.25 ml resin) into the batch incubation chamber of the spin column barrel. Wash the resin at 500 x g for 5 min.
2. Pre-equilibrate the resin with 15 ml binding buffer by centrifuging the spin

column (with the CLEAR SPIN PUSH CAP) at 750 x g for 5 min. It is critical that you repeat this step one more time with a further 15 ml fresh equilibration buffer.

NOTE: If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

CLARIFICATION OF SAMPLE

3. Pre-filter the sample through a single 0.2 µm (25 mm diameter) syringe filter.

NOTE: As with all forms of chromatography, it is critical that the sample is filtered through a final 0.2 µm syringe filter immediately before loading it on the spin column.

Optimal performance of these devices will depend on these instructions being rigorously followed.

SAMPLE LOADING

4. Transfer the spin column barrel to a fresh 50 ml centrifuge tube and load your required volume of filtered sample. The maximum sample volume is 20 ml. Tightly screw the **yellow** batch incubation cap and invert 2-3 times to mix the sample and the Ni-Advance resin. Place the column on a standard tube roller and mix for 1-3 hours.

5. After batch incubation, replace the **yellow** cap with the CLEAR spin push cap. Centrifuge the column at 750 x g for up to 10 min and collect the eluate.

NOTE: If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

Wash and Elute

6. Transfer the spin column barrel to a fresh 50 ml centrifuge tube and load your required volume of wash and elution buffer. The maximum sample volume is 20 ml.

7. Centrifuge the column at 750 x g for up to 10 min and collect the wash and eluate samples.

NOTE: If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

8. The eluate contains the target protein and is now ready for further downstream analyses.

Micro Kit Protocol

Recommended Protocol:

PRE-EQUILIBRATION

1. Pipette up to 200 μ l Ni Advance resin slurry (50% slurry) into the batch incubation chamber of the spin column barrel. Wash the resin at 13,000 x g for 20 sec.
2. Pre-equilibrate the resin with 600 μ l binding buffer by centrifuging the spin column at 13,000 x g for 20 sec. It is critical that you repeat this step one more time with a further 600 μ l fresh binding buffer.

NOTE: If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

SAMPLE LOADING

9. Transfer the spin column barrel to a fresh 2.2 ml centrifuge tube and load your required volume of filtered sample. The maximum sample volume is 600 μ l. Close the lid and vortex for 15 sec to mix the sample and resin. Repeat the vortexing every 15 min for 1 hour. Or, Place the column on a standard tube roller and mix for 15 min to 1 hour.

NOTE: If the sample needs filtering, use a Proteus Mini Clarification spin column (Protein Ark Cat # GEN-MSF500, 0.2 μ m pore size, 100 pc).

3. After batch incubation, centrifuge the column at 13,000 x g for 20 sec and collect the eluate. Spin times may need to be increased if the sample is not correctly filtered.

NOTE: If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

Wash and Elute

4. Transfer the spin column barrel to a fresh 2.2 ml centrifuge tube and load your required volume of wash and elution buffer. The maximum sample volume is 600 μ l.
5. Centrifuge the column at 13,000 x g for 20 sec and collect the wash and eluate samples.

NOTE: If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

6. The eluate contains the target protein and is now ready for further downstream analyses.

Standard IMAC buffer

Example 1: Buffers for NATIVE purifications using imidazole elution.

Lysis buffer:

50 mM Sodium phosphate buffer, 0.3 M NaCl, 10 mM imidazole, pH 8.0

Wash buffer:

50 mM Sodium phosphate buffer, 0.3 M NaCl, 20 mM imidazole, pH 8.0

Elution buffer:

50 mM Sodium phosphate buffer, 0.3 M NaCl, 500 mM imidazole, pH 8.0

Example 2: Buffers for DENATURING purifications using imidazole elution.

Denaturing lysis buffer:

50 mM Sodium phosphate buffer, 0.3 M NaCl, 8 M urea (or 6 M guanidinium hydrochloride), 10 mM imidazole, pH 8.0

Denaturing wash buffer:

50 mM Sodium phosphate buffer, 0.3 M NaCl, 8 M urea (or 6 M guanidinium hydrochloride), 20 mM imidazole, pH 8.0

Denaturing elution buffer:

50 mM Sodium phosphate buffer, 0.3 M NaCl, 8 M urea (or 6 M guanidinium hydrochloride), 500 mM imidazole, pH 8.0

Example 3: Buffers for DENATURING purifications using low pH elution.

Denaturing lysis buffer:

100 mM Sodium phosphate buffer, 10 mM tris, 8 M urea, pH 8.0

Denaturing wash buffer:

100 mM Sodium phosphate buffer, 10 mM tris, 8 M urea, pH 6.3

Denaturing elution buffer:

100 mM Sodium phosphate buffer, 10 mM tris, 8 M urea, pH 4.5

Desalting and concentrating the purified protein:

Imidazole, EDTA, or acid pH solutions do not need to be removed for your downstream application since the Fastback Ni Advance resin is resistant to EDTA (up to 20mM) and DTT (up to 20mM) for IMAC chromatography – Ni ions remain bound to ligand.

Protein refolding:

Attempts can be made to renature the target protein by dialyzing it sequentially against binding buffers containing decreasing levels of urea or by passing decreasing levels of urea in binding buffer over the washed protein bound to the Ni Advance resin and eluting the refolded protein with 300 mM imidazole (between pH 7 and pH 8). Alternatively, denatured proteins can be diluted into a large volume of buffer lacking denaturant. The dispersive effect dilutes out the denaturant resulting in the re-folding of the protein.

Questions and answers:

1. What is the shelf-life of the Fastback Ni Advance Resin?

The resin is guaranteed for 2 years after the date of manufacture provided they are stored at 2-8°C.

2. Do I need to filter the buffers prepared in my laboratory?

It is good laboratory practice to filter all buffers.

3. How should I prepare my sample for metal chelate separation?

Many chromatographic procedures demand that the sample is pre-conditioned prior to loading. We recommend that all samples are filtered to at least 0.45 µm (preferably 0.20 µm) pore size. High viscosity is mostly attributed to contaminating DNA or RNA. The intrinsic viscosity of a lysate can be reduced by either drawing it through a syringe needle several times or by adding appropriate amounts of DNase and/or RNase (5-10 µg/ml) to the lysis buffer and incubating the mix on ice for 15 mins.

4. Should I add β-mercaptoethanol (or DTT) to the lysis buffer?

Reducing agents are compatible up to 20 mM DTT and 20 mM EDTA with this resin.

5. How can I regenerate the metal chelate resin?

We recommend that you wash the resin with elution buffer and then re-equilibrate the resin with binding buffer. Proceed to the pre-equilibration step if resin is to be re-used immediately. After regeneration, the resin can also be stored in a screw-capped bottle containing 0.1 % sodium azide (made up in distilled water) at 2-8 °C until further use.

6. How can I ensure that levels of contaminants in the final eluate remain low?

We recommend that the binding buffer contains minimum 10 mM imidazole and the wash buffer contains minimum 20-30 mM imidazole.

7. Should I be concerned if the resin partially dries out during the chromatographic steps?

The resin is robust. Partially dried resin rehydrates rapidly. There are no adverse effects upon the performance of the resin.

8. Should I remove imidazole after the final elution step?

You should always remove imidazole if the protein is going to be stored. Otherwise, the protein may precipitate out of solution at -20°C or -80°C.

9. Can I load purified protein immediately on to an SDS-gel?

Proteins purified under native conditions can be loaded on to an SDS-polyacrylamide gel. Those proteins purified under denaturing conditions in 6-8 M urea can also be loaded directly on to a denaturing SDS-polyacrylamide gel.

However, proteins purified in the presence of 4-6 M guanidine hydrochloride should be buffer exchanged using an ultrafiltration device into a buffer lacking the denaturant prior to loading on to a denaturing SDS-PAGE.

10. *Do I need to remove the His-tag from the recombinant protein after purification?*

For most applications it is not necessary to remove the terminal His-tag. However, it can affect the activity, stability, or structure determination. If required a protease cleavage site (e.g. Factor Xa Protease, TEV, or enterokinase) can be engineered between the His-tag and the target protein. The tag can then be cleaved off and the protein re-purified by passing it back through the IMAC column in order to remove the digested tag and undigested His-tagged protein.

Troubleshooting assistance:

The sample does not flow easily through the spin column:

- The resin is clogged with particulates. Pre-filter the sample just before loading it on to the metal chelate resin.
- If the resin is not stored at 2-8°C, or it has been used more than once and stored in the absence of a bacteriostat, microbial growth may restrict flow through the spin column.

No elution of the target protein is observed from the resin:

- The elution conditions are too mild to desorb the target protein. Use a higher concentration of imidazole or lower the elution pH further!
- Ensure that the resin is blue in appearance. Otherwise the expressed protein will not bind effectively to the resin.
- Ensure that there are no chelators or reductants in the sample which will interfere with binding of the target protein to the resin.
- The protein may have precipitated in the column. Try denaturing conditions!
- The cell disruption method may have liberated proteolytic activities. Purify the protein under denaturing conditions if you do not need to purify an active protein.

The recovery of target protein is low:

- The His-tag maybe inaccessible. Try increasing the length of the spacer between the His-tag and protein, move the affinity tag to the other end of the protein or perform the purification under denaturing conditions.
- Ensure that the column binding capacity is proportionate to the level of expressed His-tagged protein. The target protein may pass through into the sample wash if the capacity of the column is insufficient for the level of expressed protein. If necessary add more resin to increase the capacity.
- Confirm levels of target protein by immunoassay. This will help determine if your cell disruption methods have been successful.
- The target protein may contain hydrophobic stretches which could have been toxic to the host bacterium, *E.coli*.

- Ensure that the protein is not insoluble i.e. exists in inclusion bodies and resides in the pellet. Solubilize the insoluble protein using 6-8 M urea or 4-6 M guanidine hydrochloride.
- Add further protease inhibitors to the buffers as the full-length protein may have been degraded by hydrolytic enzymes. Alternatively, reduce the time of expression, lower the temperature at which the protein is exposed or use special *E.coli* strains devoid of proteases.

The target protein elutes at an unexpected position:

- There may be an ionic interaction between the protein and metal chelate resin. You should maintain the ionic strength above 0.1 M.
- There may be hydrophobic interactions between the sample and the resin. In this instance, reduce the salt concentration and add suitable detergents or organic solvents.
- Co-purification of contaminants may occur if both the expressed protein and the contaminant have similar affinities for the matrix. In this case, a further chromatographic method such as gel filtration or ion exchange chromatography is recommended.

Glossary:

affinity chromatography - chromatographic separation based on a specific interaction between an immobilized ligand and a binding site on a macromolecule.

baculovirus – a virus vector for expression of recombinant proteins in insect cells.

chaotropic agent - a molecule which interferes with hydro-phobic interactions by disrupting the ordered structure of water molecules. Examples include urea and guanidine hydrochloride.

chelating agent – a compound such as EDTA or EGTA that is able to combine with a metal ion to form a structure with one or more rings.

cleared lysate – the soluble cell extract after the cell debris and other particulates have been removed by centrifugation.

column bed volume (CV) - the total volume occupied by the chromatographic packed bed. It is also referred to as the column volume or CV.

DL-Dithiothreitol (DTT) – reducing agent used to break disulphide bonds.

expression vector – a cloning vector intended for the foreign gene to be expressed in the host organism.

french pressure cell – a device that uses high shear forces to rupture microbial cells. The suspension is poured into a chamber, which is closed at one end by a needle valve and at the other end by a piston. Pressures of up to 16,000 lb/in² are applied by a hydraulic press against a closed needle valve. When the desired pressure is attained, the needle valve is fractionally opened to marginally relieve the pressure. The cells subsequently expand and rupture, thereby releasing the cellular components through the fractionally open valve.

freeze-thawing – a method that is sometimes used to break open cells by successive periods of slow freezing and thawing. Ice crystals are generated during the freezing stage, which

disrupt the cells when they melt during thawing. The method, however, is slow and releases a limited amount of subcellular components.

his - a 3 letter symbol for L-histidine

his-tag – a permanent affinity tag engineered into the expression vector upstream or downstream of the gene of interest to facilitate the purification of the recombinant protein. The His-tag doesn't normally have any effect upon the protein structure or function, it comprises 6-8x Histidine residues (polyhistidine) and has a molecular weight of 0.7-0.9 kDa

iminodiacetic acid (IDA) – chelating group used for immobilizing metal ions.

immobilized metal ion affinity chromatography (IMAC) – method of protein affinity purification using immobilized metal ions.

immobilized - bound to a surface, usually through covalent linkages.

inclusion bodies – quite a lot of proteins form insoluble crystalline aggregates known as inclusion bodies when they are expressed at high levels inside bacteria. The proteins can be solubilized using denaturants such as 8 M urea or 6 M guanidine hydrochloride.

ion exchange chromatography - chromatographic separation based on different charge properties of macromolecules.

isoelectric point - the pH at which the protein has no net charge.

lysozyme – an enzyme that hydrolyzes β -1,4-linkages between N-acetylmuramic acid and 2-acetamido-2-deoxy-D-glucose in peptidoglycan heteropolymers of prokaryotic cell walls. An example is egg white lysozyme and this enzyme is used to disrupt cells in order to liberate expressed proteins. 1 mg/ml lysozyme is normally added to *E.coli* cells in lysis buffer and incubated for 30 min to aid cell disruption. The pH optimum for lysozyme is pH 9.2 (Davies *et al* 1969).

nitrilotriacetic acid (NTA) – chelating group used for immobilizing metal ions.

recombinant protein – a protein coded for by a cloned gene which has often been modified to increase the expression of that protein or to alter the properties of the protein.

sonication – this technique uses ultrasonic energy to generate high transient pressures that are believed to disrupt the cells.

tris(2-chloroethyl) phosphate (TCEP) – strong irreversible reducing agent used to break disulphide bonds.

truncate - terminate prematurely or to shorten by cutting.

Literature:

Affinity Separations: A Practical Approach 1997 (Matejtschuk, P Ed.) IRL PRESS at Oxford University Press. ISBN: 0-19-963550-1.

Biochemistry, 2nd Edition 1996 (Mathews, CK & van Holde, KE) Benjamin/Cummings Publ. Co. ISBN: 0-8053-3931-0.

Cameron, I. R *et al* (1989). *Trends in Biotechnology* 7, 66-70.

Cuatrecasas, P., Wilcheck, M. and Anfinsen, C.B. (1968) *Proc. Natl. Acad. Sci. U.S.A.* 61, 636.

Davies, A., Neuberger, A. and Wilson, B. (1969). *Biochim. Biophys. Acta* 178, 294.

Protein Purification Applications: A Practical Approach 1995 (E.L.V. Harris and S. Angal Eds.) IRL PRESS at Oxford University Press. ISBN: 0-19-963023-2.

Recombinant Protein Protocols 1997 (RS Tuan Ed.) Humana Press, ISBN: 0-89603-400-3.

Porath, J., Carlson, J., Olsson, I., and Belfrage, G. (1975) *Nature* 258, 598.

Porath, J. and Olin, B. (1983) *Biochemistry* 22, 1621

Porath, J. (1988) *Trends Anal. Chem.* 7, 254.

Porath, J. (1992). *Protein Expression and Purification: Immobilized metal ion affinity chromatography* 3, 263.

Sulkowski, E. (1989) *BioEssays*, 10, 170.

Ordering information:

Product Code	Product Description
GEN-NiADV-2Micro	Proteus NiADV (Micro) 2 column kit (includes 0.2 ml NiADV resin vial)
GEN-NiADV-24Micro	Proteus NiADV (Micro) 24 column kit (includes 2.4 ml NiADV resin vial)
GEN-NiADV-100Micro	Proteus NiADV (Micro) 100 column kit (includes 10.0 ml NiADV resin vial)
GEN-NiADV-2M	Proteus NiADV (Mini) 2 column kit (includes 2 x 0.25 ml NiADV resin vials)
GEN-NiADV-12M	Proteus NiADV (Mini) 12 column kit (includes 12 x 0.25 ml NiADV resin vials)
GEN-NiADV-48M	Proteus NiADV (Mini) 48 column kit (includes 48 x 0.25 ml NiADV Resin vials)

Technical support:

Contact the Protein Ark technical support and sales centre for assistance:

Telephone +44 (0) 114 224 2236

FAX: +44 (0) 114 224 2222

Email: info@proteinark.com

Web: www.proteinark.com

Disclaimer:

- This product is for research use only and is not intended for use in clinical diagnosis. No claims beyond replacement of unacceptable material or refund of purchase price shall be allowed.