

Prime Taq DNA Polymerase

| Cat. No. | Size | Remark |
|----------|-----------|--|
| G-1000 | 250 Units | with dNTPs mix., +MgCl ₂ buffer |

Package Information

| | |
|--------|--|
| G-1000 | - Prime Taq DNA Polymerase (5 Units/μl): 50 μl - 10X Reaction buffer (with MgCl ₂): 1.0 ml - 10 mM dNTP Mixture (2.5 mM of each dNTPs): 0.5 ml |
|--------|--|

Description

Prime Taq DNA Polymerase is a high quality recombinant enzyme and catalyzes 5'→3' synthesis of DNA. The enzyme has no detectable 3'→5' proofreading exonuclease activity.

It is provided with 10X reaction buffer that contains PCR enhancers. This reaction buffer will enable or improve sub-optimal PCR caused by templates that have a high degree of secondary structure or that are GC-rich.

Buffer and Reagents

Storage Buffer

20 mM Tris-HCl (pH 8.0), 100 mM KCl, 0.5 mM EDTA, 0.1 mM DTT, 0.5% Tween 20, Nonidet P-40 and 50% Glycerol

10X Reaction Buffer (with MgCl₂)

Contains Tris-HCl (pH 9.0), 20 mM MgCl₂, (NH₄)₂SO₄ and PCR enhancers.

10 mM dNTP mixture

2.5 mM of each dATP, dCTP, dGTP and dTTP

Usage Information

- A DNA fragment which is amplified by Prime Taq DNA Polymerase has A-overhang, and it enables you to do cloning by using T-vectors.
- This product is sold for research purpose only. This is not to be used for human diagnostic or drug purposes.
- All claims must be brought within expired date.

Protocol

Optimal reaction conditions, such as reaction time, temperature and amount of template DNA, may vary and must be individually determined.

1. Thaw 10X Reaction buffer.
2. Prepare a master mix.

| Components | Volume | Final Conc. |
|-------------------------------------|--------------|--------------|
| 10X Reaction buffer | 5 μl | 1X |
| 10 mM dNTP Mix. (2.5 mM of each) | 1 ~ 5 μl | 0.2 ~ 1.0 mM |
| Upstream Primer | Variable | 0.1 ~ 1.0 μM |
| Downstream Primer | Variable | 0.1 ~ 1.0 μM |
| Prime Taq DNA Pol. (5Unit/μl) | 0.2 ~ 0.5 μl | 1.0 ~ 2.5 U |
| Template DNA | Variable | 10 fg~1 μg |
| Sterilized D.W | Variable | - |
| Total Volume | 50 μl | - |

* Amount of template:

- Bacteriophage lambda, cosmid, plasmid: 10 fg ~ 300 ng
- Total genomic DNA: 100 ng ~ 1 μg

3. Mix the master mix and dispense appropriate volumes into PCR tubes. Centrifuge the reactions in a micro-centrifuge for 10 seconds.

4. Perform PCR using your standard parameters (3-step cycling).

| Step | Temp. & Time | | Cycles |
|----------------------|-----------------|-----------------------|---------|
| | Temp. | Time | |
| Initial denaturation | 95°C | 5 min. | 1 |
| Amplification | 95°C | 30 sec. | 25 ~ 35 |
| | 50~68°C 72°C | 30 sec. 30~60 sec. | |
| Final extension | 72°C | 5 min. | 1 |

* For PCR products longer than 3~4kb, use an extension time of approximately 1 min per Kb DNA.

5. Separate the PCR products by agarose gel electrophoresis and visualize with EtBr or any other means.

■ Research Use Only

■ Store at -20°C

■ Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes.