

PrimePrep™ PCR Purification Kit

Introduction

PrimePrep™ PCR Purification Kit offer simple, rapid and cost-effective method for purification from PCR/enzyme reaction mixtures.

The purified DNA can be directly used in ligation, sequencing and other downstream application.

Kit Components

Cat. No. Reagents	K-7000 (50 prep.)	K-7001 (200 prep.)
Spin column	50 ea	50 ea x 4
Buffer PCR-B	30 ml	120 ml (40 ml x 3)
Buffer PW	10 ml	30 ml (15 ml x 2)
Buffer PE	10 ml	20 ml

Before you begin

- **Add ethanol to Buffer PW before use.**

→ Add 40 ml (K-7001: 60 ml) of absolute ethanol before use.

Experimental Protocol

- 1. Add 5 volumes of Buffer PCR-B to 1 volume of the sample and mix well by vortexing.**
If the PCR reaction product is 50 ul, add 250 ul of Buffer PCR-B.
- 2. Centrifuge the tube briefly at room temperature.**
- 3. Transfer the mixture to a Spin column.**
- 4. Centrifuge for 30 sec ~ 1 min at 13,000 rpm. Discard the flow-through and re-inserting the spin column to the collection tube.**
- 5. Add 700 ul Buffer PW and centrifuge for 30 sec. at 13,000 rpm.**
Discard the flow-through and re-inserting the spin column to the collection tube.
- 6. Centrifuge for an additional 1 ~ 2 min at 13,000 rpm to remove residual wash buffer.**
Residual ethanol of washing buffer may inhibitor subsequent enzymatic reaction.
- 7. Transfer the spin column to new 1,5 ml microcentrifuge tube.**
The 1.5 ml microcentrifuge tube is not provided.
- 8. Add 50 ul of Buffer PE or deionized distilled water to the center of the membrane in the column, let stand for 1 min and centrifuge for 1 min at 13,000 rpm.**
For larger fragment (>5kb), use pre-warmed (70°C) Buffer PE for best efficiency.