

100 bp DNA Marker

Product Name	Cat. No.	Size
100 bp DNA Marker	M-1000	0.5 ml X 1
	M-1002	0.5 ml X 5

Specification


Concentration: 103.8 ng / μ l

Recommended loading: 3.0 ~ 5.0 μ l / 5 mm lane width

Number of bands: 11

Size range: 100 ~ 1,500 bp

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA,
 5% Glycerol, 0.005% Bromophenol Blue,
 0.005% Xylene Cyanol



4.0 μ l loading / 5 mm lane width
 1.8% TBE agarose gel with stained EtBr

● Research Use Only

● Store at -20°C