

100 bp Plus DNA Marker

| Product Name | Cat. No. | Size |
|------------------------|----------|------------|
| 100 bp Plus DNA Marker | M-2000 | 0.5 ml X 1 |
| | M-2002 | 0.5 ml X 5 |

Specification


Concentration: 113 ng / μ l

Recommended loading: 4.0 ~ 5.0 μ l / 5 mm lane width

Number of bands: 18

Size range: 100 ~ 10,000 bp

Storage buffer: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA,
 5% Glycerol, 0.005% Bromophenol Blue,
 0.005% Xylene Cyanol



5.0 μ l loading / 5 mm lane width

1.0% TAE agarose gel with stained EtBr

● **Research Use Only**

● **Store at -20°C**