

# SAFETY DATA SHEET

#### **1. Identification of the Product and of the Company**

#### Identification of Product

Product Name: PwrPOPTM-7, PwrPOPTM-6, and PwrPOPTM-4 Sequencing Polymers Product Code: 037003, 037007, 037028, 037140, 037280, 037840 036005, 036010, 036003, 036007,036028 034005, 034010, 034003, 034007, 034028

#### **Company Identification**

AdvancedSeq LLC 6654 Owens Drive Pleasanton, CA 94588

Phone: +1 818-639-4362

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

# 2. Hazards Identification

# **GHS - Classification**

Signal Word WARNING



#### **Health Hazard**

Skin Corrosion/Irritation	Category 2
Serious Eye Damage/Eye Irritation	Category 2

# Physical Hazards

Not hazardous

#### **Hazard Statements**

H315 - Causes skin irritation H319 - Causes serious eye irritation

#### **Precautionary Statements**

**P305 + P351 + P338** - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing

P337 + P313 - If eye irritation persists: Get medical advice/attention

P332 + P313 - If skin irritation occurs: Get medical advice/attention P264 - Wash hands thoroughly after handling

P280 - Wear protective gloves/protective clothing/eye protection/face protection

P362 - Take off contaminated clothing and wash before reuse

Revision Date Product Code

## Principle Routes of Exposure/ Potential Health Effects

Skin:	Irritating to skin.
Inhalation:	Irritating to respiratory system.
Ingestion:	May be harmful if swallowed. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.

Specific Effects

Carcinogenic Effects: None Mutagenic Effects: None Reproductive Toxicity: None Sensitization: None Target Organ Effects: None under normal use conditions

#### HMIS

Health	1
Flammability	0
Reactivity	0

## 3. Composition/information on ingredients

Chemical Name		CAS-No EINECS-No		Weigh	
Urea		57-13-6	200-315-5	30	

We recommend handling all chemicals with caution.

# 4. First Aid Measures

Skin Contact:	Wash off immediately with plenty of water for at least 15 minutes. Remove and wash contaminated clothing before re-use. Immediate medical attention is required.
Eye Contact:	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Immediate medical attention is required.
Ingestion:	Call a physician or Poison Control Centre immediately. Never give anything by mouth to an unconscious person. Do not induce vomiting without medical advice.
Inhalation:	Move to fresh air. If not breathing, give artificial respiration. Call a physician or Poison Control Center immediately.
Notes to Physicia	n: Treat symptomatically.

#### **5. Fire-Fighting Measures**

**Suitable Extinguishing Media**: Water spray. Carbon dioxide (CO2). Foam. Dry chemical. **Special Protective Equipment for Firefighters**: Wear self-contained breathing apparatus and protective suit.

#### **6. Accidental Release Measures**

Personal Precautions	Ensure adequate ventilation. Avoid contact with skin, eyes and clothing. Use
	personal protective equipment.
Methods for Cleaning up:	Soak up with inert absorbent material.

# **Environmental Precautions**

Prevent further leakage or spillage if safe to do so. Prevent product from entering drains. See Section 12 for additional information.

## 7. Handling and Storage

**Handling:** Always wear recommended Personal Protective Equipment. No special handling advice required. **Storage:** Keep in a dry, cool and well-ventilated place. Keep in properly labelled containers.

## 8. Exposure Controls/Personal Protection

## **Exposure Limits**

Chemical Name	OSHA PEL	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Urea	none	none	none	none

Engineering measures: Ensure adequate ventilation, especially in confined areas.

#### Personal Protective Equipment

Personal Protective Equipment requirements are dependent on the user institution's risk assessment and are specific to the risk assessment for each laboratory where this material may be used.

Respiratory Protection:	In case of insufficient ventilation wear suitable respiratory equipment.
Hand Protection:	Impervious gloves.
Eye Protection:	Safety glasses with side-shields.
Skin and Body Protection:	Lightweight protective clothing.
Hygiene Measures:	Handle in accordance with good industrial hygiene and safety practice.

#### **Environmental Exposure Controls**

Prevent product from entering drains. Do not allow material to contaminate ground water system.

#### 9. Physical and Chemical Properties

#### **General Information**

Form: Liquid Appearance: No information available Odor: No information available Boiling Point/Range: No information available Melting Point/Range: No information available Flash Point: No information available Auto Ignition Temperature: No information available Oxidizing Properties: No information available. Water Solubility: soluble

#### **10. Stability and Reactivity**

Stability: Materials to Avoid: Hazardous Decomposition Products: Polymerization: Stable under normal conditions. No information available. No information available. None under normal processing.

#### **11. Toxicological Information**

#### Acute toxicity

Chemical Name Urea		Name	LD50 (oral, rat/mouse)	LD50 (dermal, rat/rabbit)	LC50 (inhalation, rat/mouse)
		l	= 8471 mg/kg (Rat)	no data available	no data available
	Revision Date Product Code	036005, 0360	Ver. 2.0 7,037028,037140,037280,037840, 110, 036003,036007, 036028, 10,034003,034007,034028	Product Name	Page 3 / 5 PwrPOP <sup>™</sup> Sequencing Polymers

#### Principle Routes of Exposure/ Potential Health Effects

Eyes: Skin: Inhalation: Ingestion:	Irritating to eyes. Irritating to skin. Irritating to respiratory system. May be harmful if swallowed Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.
Carcinogenic Effects: Mutagenic Effects: Reproductive Toxicity: Sensitization:	None None None None
Target Organ Effects:	None under normal use conditions

# **12. Ecological Information**

Ecotoxicity Effects:	No information available.
Mobility:	No information vailable.
Biodegradation:	No information vailable.
Bioaccumulation:	No information vailable.

Chemical Name	Freshwater Algae Dat	Water Flea Data	Freshwater Fish	Microtox Data	log Pow
Urea 57-13-6		Daphnia magna EC50=3910 mg/L (48 h) Daphnia magna Straus EC50>10000 mg/L			logPow-1.59

# **13. Disposal Considerations**

Dispose of in accordance with local regulations.

# **14. Transport Information**

<u>IATA</u>

Proper Shipping Name:Not classified as dangerous in the meaning of transport regulationsHazard Class:NoneSubsidiary Class:NonePacking Group:NoneUN-No:None

# **15. Regulatory Information**

Compone	TSC
Urea	Listed
57-13-6(30-60)	

# **U.S. Federal Regulations**

# **SARA 313**

This product is not regulated by SARA.

Revision Date Product Code	20-Apr-2020 Ver. 2.0 037003,03707,037028,037140,037280,037840, 036005, 036010, 036003,036007, 036028,
	034005,034010,034003,034007,034028,

# Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain HAPs.

# U.S. State Regulations

Chemical Name	Massachusetts - RTK	New Jersey - RTK	Pennsylvania - RTK	Illinois - RTK	Rhode Island - RTK
Urea	-	-	-	-	-

#### **California Proposition 65**

This product does not contain chemicals listed under Proposition 65

# **WHMIS Hazard Class**

D2B Toxic Materials



This product has been classified according to the hazard criteria of the CPR and the MSDS contains all of the information required by the CPR

# **16. Other information**

For research use only. Not intended for human or animal diagnostic or therapeutic uses. The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS MSDS DOES NOT CONSTITUTE A WARRENTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PUPOSE.

# End of Safety Data Sheet