# Protein /rk

# Datasheet

TCEP (Tris(2-carboxyethyl)phosphine hydrochloride) is an effective reducing agent for cleaving protein disulphide bonds in aqueous solution. It is soluble in water, odourless and non-toxic. It offers considerable advantages over 2-mercaptoethanol or DTT during Ni-IMAC chromatography.



### Specification:

•		
Formula:	C <sub>9</sub> H <sub>15</sub> O <sub>6</sub> P●HCl	
Molecular Weight:	286.65 g/mol	
CAS:	[51805-45-9]	
Purity:	98.0% min. (GC)	
Appearance:	White powder	
IR:	Conforms to structure	
Water content:	0.5% max. (KF)	
Solubility:	Clear to slightly hazy (50 mg/ml in water)	
Storage:	Room temperature	
Usage:	TCEP is acidic when dissolved in water with a pH of approximately 2.5. Adjust the pH of the solution with the appropriate working buffer just before use.	

# **Ordering Information:**

Product	Units	Order Code
ТСЕР	1 g	GEN-TCEP-1
ТСЕР	10 g	GEN-TCEP-10

# Protein Ark Limited

Telephone+44 (0) 114 224 2257Email:info@proteinark.com

Protein /

Calibre Scientific Unit 5A, R-Evolution, Advanced Manufacturing Park, Selden Way, Rotherham S60 5XA Email: info@proteinark.com Telephone: +4 Web: www.proteinark.com