



PRODUCT DATA SHEET

HCMV Grade 2 Antigen Solution (ready-to-use)

Catalog No.: HCMV-GR2-015

Product Name: HCMV Grade 2 Antigen Solution, inactivated by gamma radiation

Packaging Size: HCMV Grade 2 antigen dissolved in serum-free CTL-Test™ Medium, 1.5ml per vial.

Concentration: 60µg/ml

Description: This preparation contains HCMV antigens from all parts of the virus cycle, and residual host cellular material, (inactivated by gamma radiation).

Activity: HCMV stimulates antigen-specific CD4+ memory T cells.

Recommended After thawing, the reagent is a 2x ready-to-use working solution (concentration 60 µg/ml).

Usage: When working in 96 well plates, add 100µL of antigen solution to 100uL of PBMC per well. The PBMC concentration should be 100,000 – 800,000 cells/well.

Storage: Keep frozen at -80°C.

Warning: **For laboratory research use only, not for diagnostic use. Not for use in humans.**

Please note that the information contained in this document is privileged and confidential and protected from disclosure. You are hereby notified that any dissemination, distribution, or copying of this document is strictly prohibited. © 2012. All rights reserved.